



### Caring for Your SPIO

- SPIO's should be washed at least every other day to ensure proper compression.
- Machine wash in warm or cold water using mild detergent.
- Lay flat or hang to dry – DO NOT DRY IN A TUMBLE DRYER.
- When washing the TLSO be sure the front panel is attached to the back panel to prevent the Velcro from snagging the Lycra and remove X-Panels and/or Stays.
- Abuse, altering, or undue rough wear and tear will void the warranty. This includes improper washing and drying in a tumble dryer.
- The fabric is not fire resistant, so the items are not recommended for use during sleep time.

### Fitting Your SPIO

- To provide maximum benefit a physician or therapist should monitor the use and effect of the SPIO.
- When you receive the SPIO try it on first before you decide it is too small. It will look too small but may not be. If it is the same size as the child's regular t-shirt, it is too big, and you are not getting any compression.
- You should be able to place your hand between the child's body (in the trunk area) and SPIO and feel significant resistance when you pull your hand away from the child's body. An approximate one-inch resistance is adequate. More and it is too loose, less and it is too tight.
- When donning the **Upper Body Orthosis (UBO)**, unzip and if the child is able, have the child step into the UBO and then feed the arms through the sleeves. Make sure you have the sleeves turned properly so they are not twisted, as this can cause a bad fit in the armpit area. The hem of the UBO should fit over the bottom and end at the trochanter.
- When donning the **TLSO**, attach one shoulder strap and the same side panel to the neoprene back. With the child sitting or lying on their back, slide the child's arm through the fastened opening and attach the other shoulder strap. The child can be standing, sitting, or lying down on their back, depending on the child's comfort level. Place the child on their stomach on your lap or on the floor (sitting, standing or kneel-stand are also possible). To attach the side panel, pull the TLSO down over the hips and wrap the panel snugly around the child. Fasten the crotch strap so it is snug, not tight. The crotch strap is designed to prevent the TLSO from riding up. The TLSO should fit snugly and firmly. You may have to adjust it a few times after the initial donning to optimize the fit. If you slip your hand under the abdominal area, with proper tension, you should feel a snug compression. After fitting and adjustment, you may add the para-spinal stays to the back panel, if needed. To do so simply attach the stays via the Velcro to the back panel on either side of the spine, or where additional support is needed.
- When donning the **Lower Body Orthosis (LBO)**, bunch the legs like you would a pair of tights, and slide them up the child's leg. Make sure you have the seams turned properly so they are not twisted. The label or logo should be on the front.
- *After donning the SPIO, monitor the child's reaction. If they get red marks under their arms or on their skin at the sides, it may be too small, and you will want to exchange the SPIO for next larger size.*



### Caring for Your FabriFlex

- Hand or machine wash in cold water, dry in dryer on low heat (use dryer bag to reduce lint), or preferably, line dry.
- Do not use fabric softer, a brush for cleaning, or iron.
- Clean lint off fabric with a sticky lint roller after every wash.

### Using Your FabriFlex

- Using the Straps, start at the bottom of the limb and work your way up, stopping at the desired point and overlapping the fabric as you go. Add additional straps for longer length, or multiple layers.
- Using the Wraps, wrap at the desired point of support or stabilization, increasing the tension as needed.

### Warranty

Each SPIO and FabriFlex is guaranteed for fit, materials, and workmanship for 90 days. **See reverse for return/exchange information.**

**Return/Exchange Form**

**Date:** \_\_\_\_\_

1. Complete this Return/Exchange Form
2. Repackage the item and include a copy of the original Invoice or Packing Slip.
3. \*Arrange shipping to: SPIO RETURNS 127 SW 156<sup>th</sup> St. Burien WA 98166.
4. *Items must be returned with all original components for a full refund.*

***\*SPIO does not provide return postage or refund your original postage in the case of a return.***

Name: \_\_\_\_\_ Phone: \_\_\_\_\_

Address: \_\_\_\_\_

E-mail: \_\_\_\_\_

Invoice Number: \_\_\_\_\_ P.O. Number: \_\_\_\_\_

**Item(s) you are returning or exchanging:**

Item: \_\_\_\_\_ Qty: \_\_\_\_\_ Reason: \_\_\_\_\_

\_\_\_\_\_

Item: \_\_\_\_\_ Qty: \_\_\_\_\_ Reason: \_\_\_\_\_

\_\_\_\_\_

Item: \_\_\_\_\_ Qty: \_\_\_\_\_ Reason: \_\_\_\_\_

\_\_\_\_\_

**Item(s) you would like in exchange:**

Item: \_\_\_\_\_ Qty: \_\_\_\_\_

Item: \_\_\_\_\_ Qty: \_\_\_\_\_

Item: \_\_\_\_\_ Qty: \_\_\_\_\_