





Measurement Key



- Chest measurements are taken with patient lying down
- Pant measurements are taken with patient standing
- If in-between sizes refer to length measurements.
- Call for color availability

SPIO® Classic TLSO

Chest circumference is measured at nipple line.
Length is from Acromion to Greater Trochanter.

Item #	
V44	Chest Circumference: 44 cm, 30.75 cm long
V47	Chest Circumference: 47 cm, 34.5 cm long
V51	Chest Circumference: 51 cm, 38.25 cm long
V54	Chest Circumference: 54 cm, 44 cm long
V57	Chest Circumference: 57 cm, 47.75 cm long
V60	Chest Circumference: 60 cm, 51.5 cm long
V64	Chest Circumference: 64 cm, 55.25 cm long
V67	Chest Circumference: 67 cm, 59 cm long

SPIO® Lower Body Orthosis

Waist circumference is measured at belly button.
Length is from waist to Ankle (P) or bottom of Patella (SP).

Short #/Long #

SP43 / P43	Waist Circumference: 43 cm, SP 27.30 cm long / P 40.33 cm long
SP45 / P45	Waist Circumference: 45 cm, SP 30.48 cm long / P 45.00 cm long
SP47 / P47	Waist Circumference: 47 cm, SP 32.00 cm long / P 47.95 cm long
SP49 / P49	Waist Circumference: 49 cm, SP 33.02 cm / P 49.80 cm long
SP51 / P51	Waist Circumference: 51 cm, SP 34.29 cm / P 52.70 cm long
SP53 / P53	Waist Circumference: 53 cm, SP 36.83 cm / P 56.84 cm long
SP55 / P55	Waist Circumference: 55 cm, SP 39.37 cm / P 63.19 cm long
SP57 / P57	Waist Circumference: 57 cm, SP 42.54 cm / P 66.60 cm long
SP59 / P59	Waist Circumference: 59 cm, SP 45.30 cm long / P 70.04 cm long
SP61 / P61	Waist Circumference: 61 cm, SP 46.00 cm long / P 71.70 cm long
SP63 / P63	Waist Circumference: 63 cm, SP 48.89 cm long / P 77.80 cm long
SP65 / P65	Waist Circumference: 65 cm, SP 50.16 cm long / P 79.30 cm long
SP67 / P67	Waist Circumference: 67 cm, SP 54.70 cm long / P 81.70 cm long
SP69 / P69	Waist Circumference: 69 cm, SP 55.25 cm long / P 83.20 cm long

SPIO® LSO

Chest circumference is measured at nipple line.
Length is bottom of Rib Cage to Greater Trochanter.

Item #	
AB51	Chest Circumference: 51 cm, 15.5 cm long
AB54	Chest Circumference: 54 cm, 17.75 cm long
AB57	Chest Circumference: 57 cm, 21 cm long
AB60	Chest Circumference: 60 cm, 22.25 cm long
AB64	Chest Circumference: 64 cm, 25.5 cm long
AB67	Chest Circumference: 67 cm, 28 cm long
AB70	Chest Circumference: 70 cm, 32 cm long
AB73	Chest Circumference: 73 cm, 34.5 cm long
AB77	Chest Circumference: 77 cm, 37.5 cm long

SPIO® Upper Body Orthosis

Chest circumference is measured at nipple line.
Length is from Acromion to Greater Trochanter.

Short #/Long #

SS41/ S41	Chest Circumference: 41 cm, 23.00 cm long
SS44/ S44	Chest Circumference: 44 cm, 27.94 cm long
SS47/ S47	Chest Circumference: 47 cm, 31.11 cm long
SS51/ S51	Chest Circumference: 51 cm, 35.56 cm long
SS54/ S54	Chest Circumference: 54 cm, 39.00 cm long
SS57 / S57	Chest Circumference: 57 cm, 45.72 cm long
SS60/ S60	Chest Circumference: 60 cm, 47.00 cm long
SS64/ S64	Chest Circumference: 64 cm, 50.50 cm long
SS67/ S67	Chest Circumference: 67 cm, 56.40 cm long

SPIO® Arm Orthosis

Arm circumference is measured at arm pit.

Item #

AS13	Arm Circumference: 13 cm, 19.05 cm long
AS14	Arm Circumference: 14 cm, 20.06 cm long
AS16	Arm Circumference: 16 cm, 22.86 cm long
AS17	Arm Circumference: 17 cm, 24.07 cm long
AS18	Arm Circumference: 18 cm, 28.00 cm long
AS20	Arm Circumference: 20 cm, 30.48 cm long
AS21	Arm Circumference: 21 cm, 32.04 cm long
AS22	Arm Circumference: 22 cm, 35.56 cm long
AS24	Arm Circumference: 24 cm, 38.10 cm long
AS25	Arm Circumference: 25 cm, 40.00 cm long
AS26	Arm Circumference: 26 cm, 42.50 cm long
AS28	Arm Circumference: 28 cm, 45.72 cm long

SPIO® Wrist-Hand Orthosis

Palm width is measured at widest part above thumb.

Item #

G5	Width of Palm: 5 cm, 11.43 cm long
G5.5	Width of Palm: 5.5 cm, 12.19 cm long
G6	Width of Palm: 6 cm, 13.03 cm long
G6.5	Width of Palm: 6.5 cm, 14.04 cm long
G7	Width of Palm: 7 cm, 15.55 cm long
G7.5	Width of Palm: 7.5 cm, 16.50 cm long
G8	Width of Palm: 8 cm, 17.02 cm long

SPIO® Leg Orthosis

Leg circumference is measured at top of thigh.

Item #

LS19	Leg Circumference: 19 cm, 21.75 cm long
LS20	Leg Circumference: 20 cm, 24.5 cm long
LS21	Leg Circumference: 21 cm, 27.25 cm long
LS22	Leg Circumference: 22 cm, 30 cm long
LS24	Leg Circumference: 24 cm, 32.75 cm long
LS25	Leg Circumference: 25 cm, 40 cm long
LS26	Leg Circumference: 26 cm, 43.5 cm long
LS28	Leg Circumference: 28 cm, 46.75 cm long
LS29	Leg Circumference: 29 cm, 54.5 cm long
LS30	Leg Circumference: 30 cm, 56.5 cm long
LS32	Leg Circumference: 32 cm, 58.25 cm long
LS33	Leg Circumference: 33 cm, 60.25 cm long
LS34	Leg Circumference: 34 cm, 62 cm long
LS36	Leg Circumference: 36 cm, 63.75 cm long
LS37	Leg Circumference: 37 cm, 65.75 cm long