






PRODUCT GUIDE

SYSTEM	USE	FEATURES
<p>VEST</p> 	<ul style="list-style-type: none"> • Provides shoulder-trunk-hip stability and mid-line orientation • Increases body awareness • Ideal for children with poor core muscle activation, stabilization, or weakness • Excellent trial orthosis 	<ul style="list-style-type: none"> • Double layer fabric front panel • 2mm perforated, velcro sensitive detachable neoprene back panel • Adjustable hook and loop shoulder and crotch strap for anchoring • Larger sizes have removable anchor strap between the legs for toileting • Back panel may be trimmed for fit
<p>SHIRT</p>   <p>Back neck closure*</p>	<ul style="list-style-type: none"> • Assists with primary sensory awareness deficits in upper extremity and trunk • Positive impact on shoulder/arm stability and function • Improves fine motor control 	<ul style="list-style-type: none"> • Single layer of fabric • Hook and loop back neck closure • Long sleeves provide complete coverage of arms but can be cut if less coverage is needed • Made to fit snug over hips
<p>PANTS</p> 	<ul style="list-style-type: none"> • Assists with primary sensory awareness deficits in lower extremity and lower trunk • Combines with vest to improve upright stability and triple layer compression • Combines with shirt to improve full body sensory awareness and double layer compression 	<ul style="list-style-type: none"> • Single layer of fabric-lower rib to ankle • May be trimmed to desired length
<p>GLOVE</p> 	<ul style="list-style-type: none"> • Improves sensory awareness in hand and wrist • Increases stability of soft tissue structures which improves movement control 	<ul style="list-style-type: none"> • Single layer of fabric • Open fingers, fabricated thumb

For more product information, visit our website at www.spioworks.com
Email info@spioworks to inquire about custom made garments.