





# SPIO®

## Sizing Guidelines

### Measurement Key



### SPIO® Lower Body Orthosis

Waist circumference is measured at belly button.

Length is from waist to Ankle (P) or bottom of Patella (SP).

Short #/Long #

<b>SP41 / P41</b>	Waist Circumference: 41 cm, SP 25.25 cm long / P 38.25 cm long
<b>SP43 / P43</b>	Waist Circumference: 43 cm, SP 27.30 cm long / P 40.33 cm long
<b>SP45 / P45</b>	Waist Circumference: 45 cm, SP 30.48 cm long / P 45.00 cm long
<b>SP47 / P47</b>	Waist Circumference: 47 cm, SP 32.00 cm long / P 47.95 cm long
<b>SP49 / P49</b>	Waist Circumference: 49 cm, SP 33.02 cm / P 49.80 cm long
<b>SP51 / P51</b>	Waist Circumference: 51 cm, SP 34.29 cm / P 52.70 cm long
<b>SP53 / P53</b>	Waist Circumference: 53 cm, SP 36.83 cm / P 56.84 cm long
<b>SP55 / P55</b>	Waist Circumference: 55 cm, SP 39.37 cm / P 63.19 cm long
<b>SP57 / P57</b>	Waist Circumference: 57 cm, SP 42.54 cm / P 66.60 cm long
<b>SP59 / P59</b>	Waist Circumference: 59 cm, SP 45.30 cm long / P 70.04 cm long
<b>SP61 / P61</b>	Waist Circumference: 61 cm, SP 46.00 cm long / P 71.70 cm long
<b>SP63 / P63</b>	Waist Circumference: 63 cm, SP 48.89 cm long / P 77.80 cm long
<b>SP65 / P65</b>	Waist Circumference: 65 cm, SP 50.16 cm long / P 79.30 cm long
<b>SP67 / P67</b>	Waist Circumference: 67 cm, SP 54.70 cm long / P 81.70 cm long
<b>SP69 / P69</b>	Waist Circumference: 69 cm, SP 55.25 cm long / P 83.20 cm long
<b>SP71 / P71</b>	Waist Circumference: 71 cm, SP 57.50 cm long / P 86.20 cm long
<b>SP73 / P73</b>	Waist Circumference: 73 cm, SP 59.50 cm long / P 90.50 cm long

- Chest measurements are taken with patient lying down
- Pant measurements are taken with patient standing
- If in-between sizes refer to length measurements.
- Call for color availability

### SPIO® Classic TLSO

Chest circumference is measured at nipple line.  
Length is from Acromion to Greater Trochanter.

Item #

<b>V41</b>	Chest Circumference: 41 cm, 27 cm long
<b>V44</b>	Chest Circumference: 44 cm, 30.75 cm long
<b>V47</b>	Chest Circumference: 47 cm, 34.5 cm long
<b>V51</b>	Chest Circumference: 51 cm, 38.25 cm long
<b>V54</b>	Chest Circumference: 54 cm, 44 cm long
<b>V57</b>	Chest Circumference: 57 cm, 47.75 cm long
<b>V60</b>	Chest Circumference: 60 cm, 51.5 cm long
<b>V64</b>	Chest Circumference: 64 cm, 55.25 cm long
<b>V67</b>	Chest Circumference: 67 cm, 59 cm long
<b>V70</b>	Chest Circumference: 70 cm, 62.75 cm long
<b>V73</b>	Chest Circumference: 73 cm, 66.5 cm long
<b>V77</b>	Chest Circumference: 77 cm, 70.25 cm long

### SPIO® LSO

Chest circumference is measured at nipple line.  
Length is bottom of Rib Cage to Greater Trochanter.

Item #

<b>AB51</b>	Chest Circumference: 51 cm, 15.5 cm long
<b>AB54</b>	Chest Circumference: 54 cm, 17.75 cm long
<b>AB57</b>	Chest Circumference: 57 cm, 21 cm long
<b>AB60</b>	Chest Circumference: 60 cm, 22.25 cm long
<b>AB64</b>	Chest Circumference: 64 cm, 25.5 cm long
<b>AB67</b>	Chest Circumference: 67 cm, 28 cm long
<b>AB70</b>	Chest Circumference: 70 cm, 32 cm long
<b>AB73</b>	Chest Circumference: 73 cm, 34.5 cm long
<b>AB77</b>	Chest Circumference: 77 cm, 37.5 cm long

### SPIO® Upper Body Orthosis

Chest circumference is measured at nipple line.  
Length is from Acromion to Greater Trochanter.

Short #/Long #

<b>SS41/ S41</b>	Chest Circumference: 41 cm, 23.00 cm long
<b>SS44/ S44</b>	Chest Circumference: 44 cm, 27.94 cm long
<b>SS47/ S47</b>	Chest Circumference: 47 cm, 31.11 cm long
<b>SS51/ S51</b>	Chest Circumference: 51 cm, 35.56 cm long
<b>SS54/ S54</b>	Chest Circumference: 54 cm, 39.00 cm long
<b>SS57 / S57</b>	Chest Circumference: 57 cm, 45.72 cm long
<b>SS60/ S60</b>	Chest Circumference: 60 cm, 47.00 cm long
<b>SS64/ S64</b>	Chest Circumference: 64 cm, 50.50 cm long
<b>SS67/ S67</b>	Chest Circumference: 67 cm, 56.40 cm long
<b>SS70/ S70</b>	Chest Circumference: 70 cm, 60.50 cm long
<b>SS73/ S73</b>	Chest Circumference: 73 cm, 64.30 cm long
<b>SS77/ S77</b>	Chest Circumference: 77 cm, 68.50 cm long

### SPIO® Arm Orthosis

Arm circumference is measured at arm pit.

Item #

<b>AS13</b>	Arm Circumference: 13 cm, 19.05 cm long
<b>AS14</b>	Arm Circumference: 14 cm, 20.06 cm long
<b>AS16</b>	Arm Circumference: 16 cm, 22.86 cm long
<b>AS17</b>	Arm Circumference: 17 cm, 24.07 cm long
<b>AS18</b>	Arm Circumference: 18 cm, 28.00 cm long
<b>AS20</b>	Arm Circumference: 20 cm, 30.48 cm long
<b>AS21</b>	Arm Circumference: 21 cm, 32.04 cm long
<b>AS22</b>	Arm Circumference: 22 cm, 35.56 cm long
<b>AS24</b>	Arm Circumference: 24 cm, 38.10 cm long
<b>AS25</b>	Arm Circumference: 25 cm, 40.00 cm long
<b>AS26</b>	Arm Circumference: 26 cm, 42.50 cm long
<b>AS28</b>	Arm Circumference: 28 cm, 45.72 cm long

### SPIO® Wrist-Hand Orthosis

Palm width is measured at widest part above thumb.

Item #

<b>G5</b>	Width of Palm: 5 cm, 11.43 cm long
<b>G5.5</b>	Width of Palm: 5.5 cm, 12.19 cm long
<b>G6</b>	Width of Palm: 6 cm, 13.03 cm long
<b>G6.5</b>	Width of Palm: 6.5 cm, 14.04 cm long
<b>G7</b>	Width of Palm: 7 cm, 15.55 cm long
<b>G7.5</b>	Width of Palm: 7.5 cm, 16.50 cm long
<b>G8</b>	Width of Palm: 8 cm, 17.02 cm long

### SPIO® Leg Orthosis

Leg circumference is measured at top of thigh.

Item #

<b>LS19</b>	Leg Circumference: 19 cm, 21.75 cm long
<b>LS20</b>	Leg Circumference: 20 cm, 24.5 cm long
<b>LS21</b>	Leg Circumference: 21 cm, 27.25 long
<b>LS22</b>	Leg Circumference: 22 cm, 30 cm long
<b>LS24</b>	Leg Circumference: 24 cm, 32.75 cm long
<b>LS25</b>	Leg Circumference: 25 cm, 40 cm long
<b>LS26</b>	Leg Circumference: 26 cm, 43.5 cm long
<b>LS28</b>	Leg Circumference: 28 cm, 46.75 cm long
<b>LS29</b>	Leg Circumference: 29 cm, 54.5 cm long
<b>LS30</b>	Leg Circumference: 30 cm, 56.5 cm long
<b>LS32</b>	Leg Circumference: 32 cm, 58.25 cm long
<b>LS33</b>	Leg Circumference: 33 cm, 60.25 cm long
<b>LS34</b>	Leg Circumference: 34 cm, 62 cm long
<b>LS36</b>	Leg Circumference: 36 cm, 63.75 cm long
<b>LS37</b>	Leg Circumference: 37 cm, 65.75 cm long