



Measurement Key



- Chest measurements are taken with patient lying down
- If in-between sizes refer to length measurements.
- Available in Black

QUEST

TLSO with Rigid Stays

Chest circumference is measured at nipple line.
Length is from Acromion to Greater Trochanter.

Item #	
VS44	Chest Circumference: 44 cm, 30.75 cm long
VS47	Chest Circumference: 47 cm, 34.5 cm long
VS51	Chest Circumference: 51 cm, 38.25 cm long
VS54	Chest Circumference: 54 cm, 44 cm long
VS57	Chest Circumference: 57 cm, 47.75 cm long
VS60	Chest Circumference: 60 cm, 51.5 cm long
VS64	Chest Circumference: 64 cm, 55.25 cm long
VS67	Chest Circumference: 67 cm, 59 cm long

*** THIS IS NOT A CUSTOM MEASUREMENT FORM ***

The CUSTOM FORMS can be found at
spioworks.com/pages/customspio

ODYSSEY

TLSO with Thermoplastic Back Panel

Chest circumference is measured at nipple line.
Length is from Acromion to Greater Trochanter.

Item #	
VT44	Chest Circumference: 44 cm, 30.75 cm long
VT47	Chest Circumference: 47 cm, 34.5 cm long
VT51	Chest Circumference: 51 cm, 38.25 cm long
VT54	Chest Circumference: 54 cm, 44 cm long
VT57	Chest Circumference: 57 cm, 47.75 cm long
VT60	Chest Circumference: 60 cm, 51.5 cm long
VT64	Chest Circumference: 64 cm, 55.25 cm long
VT67	Chest Circumference: 67 cm, 59 cm long

EXPEDITION

TLSO with X Panel and Rigid Stays

Chest circumference is measured at nipple line.
Length is from Acromion to Greater Trochanter.

Item #	
VX44	Chest Circumference: 44 cm, 30.75 cm long
VX47	Chest Circumference: 47 cm, 34.5 cm long
VX51	Chest Circumference: 51 cm, 38.25 cm long
VX54	Chest Circumference: 54 cm, 44 cm long
VX57	Chest Circumference: 57 cm, 47.75 cm long
VX60	Chest Circumference: 60 cm, 51.5 cm long
VX64	Chest Circumference: 64 cm, 55.25 cm long
VX67	Chest Circumference: 67 cm, 59 cm long