





### Measurement Key



- Chest measurements are taken with patient lying down
- Pant measurements are taken with patient standing
- If in-between sizes refer to length measurements.
- Call for color availability

### SPIO® Classic TLSO

Chest circumference is measured at nipple line.  
Length is from Acromion to Greater Trochanter.

Item #	
V44	Chest Circumference: 44 cm, 30.75 cm long
V47	Chest Circumference: 47 cm, 34.5 cm long
V51	Chest Circumference: 51 cm, 38.25 cm long
V54	Chest Circumference: 54 cm, 44 cm long
V57	Chest Circumference: 57 cm, 47.75 cm long
V60	Chest Circumference: 60 cm, 51.5 cm long
V64	Chest Circumference: 64 cm, 55.25 cm long
V67	Chest Circumference: 67 cm, 59 cm long

**\* THIS IS NOT A CUSTOM MEASUREMENT FORM \***

The CUSTOM FORMS can be found at  
[spioworks.com/pages/customspio](http://spioworks.com/pages/customspio)

### SPIO® Lower Body Orthosis

Waist circumference is measured at belly button.  
Length is from waist to Ankle (P) or bottom of Patella (SP).

Short #/Long #

SP43 / P43	Waist Circumference: 43 cm, SP 27.30 cm long / P 40.33 cm long
SP45 / P45	Waist Circumference: 45 cm, SP 30.48 cm long / P45.00 cm long
SP47 / P47	Waist Circumference: 47 cm, SP 32.00 cm long / P 47.95 cm long
SP49 / P49	Waist Circumference: 49 cm, SP 33.02 cm / P 49.80 cm long
SP51 / P51	Waist Circumference: 51 cm, SP 34.29 cm / P 52.70 cm long
SP53 / P53	Waist Circumference: 53 cm, SP 36.83 cm / P 56.84 cm long
SP55 / P55	Waist Circumference: 55 cm, SP 39.37 cm / P 63.19 cm long
SP57 / P57	Waist Circumference: 57 cm, SP 42.54 cm / P 66.60 cm long
SP59 / P59	Waist Circumference: 59 cm, SP 45.30 cm long / P 70.04 cm long
SP61 / P61	Waist Circumference: 61 cm, SP 46.00 cm long / P 71.70 cm long
SP63 / P63	Waist Circumference: 63 cm, SP 48.89 cm long / P 77.80 cm long
SP65 / P65	Waist Circumference: 65 cm, SP 50.16 cm long / P 79.30 cm long
SP67 / P67	Waist Circumference: 67 cm, SP 54.70 cm long / P 81.70 cm long

### SPIO® Upper Body Orthosis

Chest circumference is measured at nipple line.  
Length is from Acromion to Greater Trochanter.

Short #/Long #

SS44/ S44	Chest Circumference: 44 cm, 27.94 cm long
SS47/ S47	Chest Circumference: 47 cm, 31.11 cm long
SS51/ S51	Chest Circumference: 51 cm, 35.56 cm long
SS54/ S54	Chest Circumference: 54 cm, 39.00 cm long
SS57/ S57	Chest Circumference: 57 cm, 45.72 cm long
SS60/ S60	Chest Circumference: 60 cm, 47.00 cm long
SS64/ S64	Chest Circumference: 64 cm, 50.50 cm long
SS67/ S67	Chest Circumference: 67 cm, 56.40 cm long

### SPIO® LSO

Chest circumference is measured at nipple line.  
Length is bottom of Rib Cage to Greater Trochanter.

Item #

<b>AB57</b>	Chest Circumference: 57 cm, 21 cm long
<b>AB60</b>	Chest Circumference: 60 cm, 22.25 cm long
<b>AB64</b>	Chest Circumference: 64 cm, 25.5 cm long
<b>AB67</b>	Chest Circumference: 67 cm, 28 cm long
<b>AB70</b>	Chest Circumference: 70 cm, 32 cm long
<b>AB73</b>	Chest Circumference: 73 cm, 34.5 cm long
<b>AB77</b>	Chest Circumference: 77 cm, 37.5 cm long

### SPIO® Wrist-Hand Orthosis

Palm width is measured at widest part above thumb.

Item #

<b>G5</b>	Width of Palm: 5 cm, 11.43 cm long
<b>G5.5</b>	Width of Palm: 5.5 cm, 12.19 cm long
<b>G6</b>	Width of Palm: 6 cm, 13.03 cm long
<b>G6.5</b>	Width of Palm: 6.5 cm, 14.04 cm long
<b>G7</b>	Width of Palm: 7 cm, 15.55 cm long
<b>G7.5</b>	Width of Palm: 7.5 cm, 16.50 cm long
<b>G8</b>	Width of Palm: 8 cm, 17.02 cm long

### SPIO® Arm Orthosis

Arm circumference is measured at arm pit.

Item #

<b>AS13</b>	Arm Circumference: 13 cm, 19.05 cm long
<b>AS14</b>	Arm Circumference: 14 cm, 20.06 cm long
<b>AS16</b>	Arm Circumference: 16 cm, 22.86 cm long
<b>AS17</b>	Arm Circumference: 17 cm, 24.07 cm long
<b>AS18</b>	Arm Circumference: 18 cm, 28.00 cm long
<b>AS20</b>	Arm Circumference: 20 cm, 30.48 cm long
<b>AS21</b>	Arm Circumference: 21 cm, 32.04 cm long
<b>AS22</b>	Arm Circumference: 22 cm, 35.56 cm long
<b>AS24</b>	Arm Circumference: 24 cm, 38.10 cm long
<b>AS25</b>	Arm Circumference: 25 cm, 40.00 cm long
<b>AS26</b>	Arm Circumference: 26 cm, 42.50 cm long

### FabriFlex Wrap

Wraps can be used for many applications. Smaller wraps are good for younger children and/or upper limbs. Larger wraps are good for older children and/or lower limbs.

Item #

<b>F720</b>	7 inches wide x 20 inches long
<b>F930</b>	9 inches wide x 30 inches long
<b>F1135</b>	11 inches wide x 35 inches long

### FabriFlex Straps

Straps can be used for many applications. Smaller straps are good for younger children and/or upper limbs. Larger straps are good for older children and/or lower limbs.

Item #

<b>FP2</b>	Two inch 3 pack (2" wide x 30", 40", & 50" long)
<b>FP4</b>	Four inch 2 pack (4" wide x 30" & 50" long)