



SPIO®

Upper Body Orthosis Custom Measurement Form

Date _____

Name (Patient) _____ Age _____

Contact Name _____


Contact Phone _____ Contact Email _____

Measurement Key

Length = Circumference =

Instructions: Measure child laying down. Measure elbow circumference at 15 degrees flexion. Measure in centimeters. Measure both arms. Measure over diaper or underwear. Measure on top of the torso, not on the side.

ALL MEASUREMENT FIELDS ARE REQUIRED



Shoulder Tip to Shoulder Tip

 cm

Line C - D

Shoulder to Chest (nipple line)

 cm

Line E - F

Shoulder to Waist

 cm

Line E - G

Shoulder to Trochanter

 cm

Line E - H

Upper Arm to Elbow Bend

 cm

Line I - J

Elbow Bend to Wrist

 cm

Line J - K

Product Measurements

 cmArmpit to Shirt End
Line F - H

Sleeve Inseam Length

 cm

Long Sleeve

Short Sleeve

Neck Closure

Zipper

Velcro

None

Quantity

Additional Information

Abdominal Reinforcement

Yes

No

Please include a copy of the SPIO Order Form, or your purchase order, along with your custom measurement form.